

M7

# Product Launch Meeting

Product Overview



**mindray**


# MINDRAY M7

**mindray**

***Clarity . Accuracy . Speed***

***Bring Your Crystal Ball to Point-of-Care***

The Mindray M7 Diagnostic Ultrasound System is a powerful imaging tool to assist you in meeting your clinical challenges today and tomorrow. The M7 is designed for use in all point of care environments. It delivers premium imaging performance across a broad range of specialties. Providing accurate data with speed, the M7 enables clinicians to achieve enhanced level of diagnostic confidence and efficiency.



# MINDRAY M7

**mindray**

- **System Performance Improvements**
- Image Quality Advancements
- Intelligent Workflow Improvements
- Data Management Improvements
- Mobile Trolley



# System Performance Improvements **mindray**

- **175%** increase in processing power versus M5
- Achievable through MINDRAY Proprietary Chip Technologies
- System Display
  - SHARP Next Generation LCD
  - Anti-glare coating
  - **25%** brighter than M5 display
  - Improves contrast resolution
- New power design and improve battery life **150%**



# System Specification

**mindray**

- **Weight:** 14.3 lb (6.5 kg)
- **Dimensions:** 14.06" L x 14.21" W x 2.95" H  
(35.7 cm L x 36.1 cm W x 7.5 cm H)
- **Display:** 15" Diagonal LCD (NTSC or PAL)
- **Platform:** Windows XP Professional SP3
- **CPU:** INTEL Core 2 Duo
- **MEMORY:** 2GB DDR3 RAM
- **Hard Disk:** 160G SATAII HHD
- **Power:** System operates via rechargeable lithium-ion battery or AC power (universal power adapter, 100-240 VAC, 50/60 Hz input, 15 VDC output) , 1.5 - 3 hours battery life



# System Specification

**mindray**

➤ **Multi-Language Support:** English, French, German, Italian, Portuguese, Russian, Spanish, Polish, Czech, Turkish, Finnish, Danish, Icelandic, Norwegian, Swedish and Simplified Chinese. (Display and Input)



# MINDRAY M7

**mindray**

- System Performance Improvements
- **Image Quality Advancements**
- Intelligent Workflow Improvements
- Data Management Improvements
- Mobile Trolley



# Imaging Modes

**mindray**

- B
- M
- Free Xros M (Anatomical M)
- Color M (CM)
- Color Doppler
- Power (DirPower)
- Pulsed Wave Doppler
- CW Doppler
- Smart 3D
- Static 3D
- 4D
- iScape (panoramic imaging)
- TDI (Tissue Doppler imaging)



# Image Quality Technologies



-HPRF

-THI

-Phase Inversion Tissue Harmonic Imaging

-iBeam™ Spacial compound imaging

-iClear™ Adaptive Speckle Suppression Imaging

-iTouch™ Intelligent Image Optimization

-iZoom™ Automatically expand the image to full screen






# Image Quality Advancements



## Phase Inversion Tissue Harmonic Imaging

Patent technology to extract the nonlinear harmonic component without interference of fundamental wave

- Increase noise reject capability
- Better contrast resolution
- Better visualize tissue subtleties
- Support Linear, Convex and Phase Array Transducers





# Image Quality Advancements







## iBeam™ Spacial compound imaging

- New algorithm for compound imaging
- Support Linear, Convex Transducers

Before

After




# Image Quality Advancements



## iClear™ Adaptive Speckle Suppression Imaging

- Reduce Speckle Noise
- Enhanced margin definition and tissue differentiation
- Improved Contrast Resolution




# Image Quality Advancements



## iTouch™ Intelligent Image Optimization

- B/THI : Gain, TGC, Dynamic Range
- PW/CW : Baseline, PRF
- Auto Focus Adjustment



# Image Quality Advancements



- Significant image quality improvements
- Up to **50%** increase in image quality vs. M5

## Emission

- Physical Channels Upgrade **100%** vs. M5
- Wave increasing **167%** reliability

## Receiving


- Physical Channels Upgrade **100%** vs. M5
- SNR upgrade **20%** vs. Competitors

## Transducer

- Support **192** element transducer
- Frequency 2Mhz~15Mhz



# M7 Abdominal Performance

**mindray**



# M7 OB/GYN Performance



**mindray**








# M7 OB/GYN Performance

**mindray**




# M7 Cardiac Performance




# M7 Cardiac Performance

**mindray**



# M7 Vascular Performance



# M7 SM Performance




The image displays three overlapping ultrasound B-mode frames from a Mindray M7 SM system. Each frame shows a different exam type and includes technical parameters and a color scale.

- Top Frame (Thyroid):**
  - Header: Mindray Medical, 25/01/2010 17:26:14, AP 100% MI 0.6 TIS 0.0
  - Case ID: 20100125-172539-6F22
  - Probe: 7L4s
  - Exam: Thyroid
  - Mode: B1
  - Parameters: F10 / D3.7, G50 / FR24, IP5 / DR85
  - Color Scale: 11.4 to -11.4
- Middle Frame (Thyroid):**
  - Header: Mindray Medical, 25/01/2010 17:27:57, AP 100% MI 0.4 TIS 0.1
  - Case ID: 20100125-172539-6F22
  - Probe: 7L4s
  - Exam: Thyroid
  - Mode: B2
  - Parameters: F10 / D3.7, G54 / FR21, IP5 / DR85
  - Mode C: F5.7 / G57, IP5 / WF277, PRF1.7k
  - Color Scale: 11.4 to -11.4
- Bottom Frame (Breast):**
  - Header: Mindray Medical, 26/10/2009 14:47:27, AP 100% MI 0.7 TIS 0.1
  - Case ID: 20091026-143951-E876
  - Probe: 7L4s
  - Exam: Breast
  - Mode: B1
  - Parameters: F10 / D3.7, G62 / FR93, IP5 / DR55
  - Scale: 1, 2, 3
  - Page: 282/282







# M7 Anesthesia Performance









# Transducers

Transducer	Applications	Bandwidth	Depth	Biopsy Kit
 <p>C5-2s</p>	<p>Abdomen, Gynecology and Obstetrics, vascular, musculoskeletal, pediatrics, regional nerve block, procedure and biopsy guidance.</p>	<p>2-5MHz R50mm broadband Convex array</p>	<p>38.8cm</p>	<p>Available NGB-015</p>
 <p>7L4s</p>	<p>Breast, Small parts, musculoskeletal, pediatrics, abdomen, vascular, venous access, nerve block, musculoskeletal, superficial and IMT</p>	<p>4-10MHz L38mm broadband Linear array</p>	<p>28.6cm</p>	<p>Available NGB-007</p>
 <p>L14-6s</p>	<p>Superficial, small parts, vascular, musculoskeletal, pediatrics, Nerve, venous access and biopsy guidance</p>	<p>6-14MHz 25mm broadband Linear array</p>	<p>28.6 cm</p>	<p>Available NGB-016</p>
 <p>P4-2s</p>	<p>Cardiac, abdominal, obstetrics, transcranial, pediatrics</p>	<p>2-4MHz 23mm broadband phased array</p>	<p>28.6cm</p>	<p>Available NGB-011</p>

# Transducers



Transducer	Applications	Bandwidth	Depth	Biopsy Kit
 <p>P7-3s</p>	<p>Cardiac, abdominal, obstetrics, transcranial, pediatrics</p>	<p>3-7MHz 21 mm broadband phased array</p>	<p>28.6cm</p>	<p>None</p>
 <p>V10-4s</p>	<p>Obstetrics, gynecology, urology, procedure and biopsy guidance</p>	<p>4-10MHz R10mm broadband Convex array</p>	<p>28.6 cm</p>	<p>Available NGB-004</p>
 <p>V10-4Bs</p>	<p>Obstetrics, gynecology, urology, procedure and biopsy guidance</p>	<p>4-10MHz R10mm broadband Convex array</p>	<p>28.6 cm</p>	<p>Available NGB-004</p>
 <p>4CD4s</p>	<p>Abdomen, gynecology, obstetrics</p>	<p>2-6MHz R40 mm broadband Convex array</p>	<p>28.6cm</p>	<p>None</p>

# Application Specific Calculations **mindray**

Dedicated Software Package

Abdomen, Obstetrics, Gynecology, Cardiology, Small parts, Urology, Peripheral Vascular, Pediatrics, Nerve Blocks, and Emergency package



# MINDRAY M7

**mindray**

- System Performance Improvements
- Image Quality Advancements
- **Intelligent Workflow Improvements**
- Data Management Improvements
- Mobile Trolley and Peripherals



# Intelligent Workflow Improvements **mindray**

## New Control Panel Layout

- Label changes to User-Defined keys
  - improves workflow and ease of use
- New trackball design
  - Type: Capacitive
  - reduces yellowing and glare
- New White LED backlighting
  - Improves visibility by 25%





# Intelligent Workflow Improvements **mindray**





# Intelligent Workflow Improvements **mindray**



**Patient Information**

**Optimize Workflow for Vascular**

**Flexible Design**

**Volume Adjust**

**Flexible Design**

# Intelligent Workflow Improvements **mindray**

## New Widget

- Intelligent Cursor



A screenshot of a patient information form titled "General Information". The form contains several input fields: "Patient ID" with the value "090616-011016-E7EC"; "Last Name" with "ZHANGSAN"; "First Name" (empty); "M.I." (empty); "DOB" with a date picker showing "DD/MM/YYYY" and the number "9"; "Age" with a numeric input field containing "1" and a "Years" dropdown menu; and "Gender" with a dropdown menu showing "Unknown". A green rectangular cursor is positioned over the "Age" input field.



- Scroll Bar for ease of use



# Intelligent Workflow Improvements **mindray**


## Comments (Annotations)

- Easier, Faster!
  - Press the <Comment> key, and the cursor becomes “|”.
  - Press any alphanumeric key, and the corresponding letter or numeral is displayed besides the cursor.
  - **Move the Cursor, End of Editing.**



# Intelligent Workflow Improvements **mindray**

## Comments (Annotations)


- Home Position
  - “Set Home” for commands cursor positioning.
  - Auto Complete from Text Library Function
  - **Move the Cursor, End of Editing.**
  
- Add Arrow 
  - Press the <Arrow> key, Adjust the arrow, Press <Set> to anchor the arrow position.
  - Repeat the above steps to add more arrows if necessary.



# Intelligent Workflow Improvements **mindray**

## Body Marks (Pictograms)

- The system can be configured with body mark libraries (**128** total) including Abdomen, Cardiology, GYN (Gynecology), OB (Obstetrics), Urology, SMP (Small Part), Vascular, EM (Emergency) and Nerve.
- Users can configure body marks for specific exam modes, **Customize** the body marks and **import/ export** the body marks.)



# Intelligent Workflow Improvements **mindray**

## iZoom™ Automatically expand the image to full screen

- Full screen display: 1024 x 768
- One Button to get larger image, bigger fonts and increase ease of use.
- According to the region to be zoomed, the system supports two types of full-screen zooming:



# MINDRAY M7


**mindray**

- System Performance Improvements
- Image Quality Advancements
- Intelligent Workflow Improvements
- **Data Management Improvements**
- Mobile Trolley



## Cine Review Enhancement

- Cine Compare
  - The system supports to make a comparison of different frames in one cine file, and you can perform measurement, comments adding and body marks on the frame.
  - Tips: cine compare can only be performed on single-format 2D images only



## Live Capture

- Live capture refers to save images or cines in image scanning status, and after the storage, the system continues image scanning.
- Types: retrospective and prospective saving.



The screenshot shows a 'Live Capture' settings menu with two sections: 'Prospective Cine Length' and 'Retrospective Cine Length'. Each section has two rows of settings: 'Time(s)' and 'Beat(heart cycle)'. The values are displayed in a slider control with a range in parentheses.

Section	Parameter	Value	Range
Prospective Cine Length	Time(s)	180	(1-480s)
	Beat(heart cycle)	120	(1-120 Cycles)
Retrospective Cine Length	Time(s)	5	(1-120s)
	Beat(heart cycle)	120	(1-120 Cycles)

## Task Management

- iRoam™, 802.11b/g wireless data transfer solution
  - USB
    - PC format transfer: JPG/ AVI, BMP/ AVI, TIFF/ AVI.
    - DCM format transfer
    - Export report in RTF format.
  - Ethernet
    - Network Storage
    - Network Printer
  - iDock (IOM-21)
    - USB, ECG, Serial, Audio Out L/R, Mic In, Remote control, Composite Video out, DVI-I Out



# Data Management Improvements



## Connectivity Enhancement

### ➤ Support Backstage Task

➤ The system supports three types of task management:

- Storage Task: displays the DICOM storage task.
- Print Task: displays the DICOM print task
- Media Storage Task
  - DICOM media storage task
  - Backup task (system-relevant format)
  - Send to external devices
  - Network Storage task.





## DICOM 3.0 Compatible

- DICOM Basic
  - DICOM connection verify, DICOM task management, DICOM storage, DICOM print, DICOM storage commitment, DICOM media storage (DICOMDIR review)
- DICOM Worklist
- DICOM Query/Retrieve
- DICOM MPPS (Modality Performed Procedure Step)
- Structured Report (SR)
  - DICOM OB/GYN S/R
  - DICOM VAS S/R
  - DICOM CARDIAC S/R



*Digital Imaging and Communications in Medicine*

# MINDRAY M7

**mindray**


- System Performance Improvements
- Image Quality Advancements
- Intelligent Workflow Improvements
- Data Management Improvements
- **Mobile Trolley**



# Mobile Trolley

**mindray**

- Patent Design
- All-In-One Design
- Ergonomic and ease of Use
- Customer Oriented



# Mobile Trolley

**mindray**



Patent Design

---

Open/Close Drawer Design

---

Quick Lock / Unlock Design

---



# Mobile Trolley

**mindray**

Transducer Holder x 6



Probe extend module



Compartment for placing

1. ECG module and B/W printer
2. Color Printer
3. DVD RW





# Mobile Trolley

**mindray**



**Thank you!**

**mindray**