

HyPort P20/30

Medical Supply Unit

Streamlined Design
Effortless Accessibility



HyPort P20/30

Transforming clutter into order, suitable for clinical settings including operating rooms, induction rooms, ICUs, and research laboratories.

Saving space while reducing interference between devices, making veterinarian work more efficient.



HyPort P20



HyPort P30

Safety Assurance

Strictly designed following ISO 11197:2019 medical pendent standard

Electrical safety



Undergoing rigorous testing for resistance, current leakage, and grounding impedance following IEC 60601-1 standards for medical equipment.



Unique bottom exhaust structure design to meet the requirements of the leakage test in the standard *.

Fire protection

Adheres strictly to the flame retardant standard of UL 94, with a shell fire prevention rating of V-0.

Infection control

Antimicrobial powder coating surface treatment with SGS certificate. Invisible screw design for easier cleaning.

*With an oxygen flow rate of 1L/min leaking into the chamber, the oxygen concentration tested after 10 minutes should not exceed 25%.

Efficient Clinical Use

The HyPort series has a wide range of accessories that are tailored to the needs of the clinic, solving the problems of storage and cable management, and bringing a more efficient and cleaner operating room environment.

Smart storage

○ Rail-drawer with dividing grids



○ Basket



○ Endoscopy camera holder



○ Rail-mount manager box



Cable management

○ Power cable manager for pumps



○ Power cable manager



○ Leadwire manager



Highly expandable rail and side rail design

With four internal guide rails, the accessories can be mounted on both sides of the column, which could cope with complex clinical scenes.



Model	Configuration			
	Arm		Type	Column
HyPort P30 Anesthesia Pendent	Single Arm	750mm	Anesthesia	800mm
	Double Arms	750+750mm	Anesthesia	800mm
HyPort P20 Anesthesia Pendent	Single Arm	750mm	Anesthesia	800mm
	Double Arms	750+750mm	Anesthesia	600mm
HyPort P20 Surgical Pendent	Single Arm	750mm	Surgical	800mm
	Double Arms	750+750mm	Surgical	600mm



LinkedIn



Facebook



Instagram



| Search Mindray Animal Medical

www.mindrayanimal.com
service@mindrayanimal.com

P/N:ENG-HyPort P20/30-210285X4P-20240613

©2022 Shenzhen Mindray Animal Medical Technology Co., LTD. All rights reserved.

