



OP-System X-Top

Veterinary Treatment and Surgery Table

- High-End veterinary special surgery table for ophthalmology and advanced orthopedics
- Lifting column chassis of highest stability
- Individually removable and re-installable laparoscopy cut-outs on both sides (also allows for closer access to the patient)
- Developed for complex procedures and surgery - also possible while seating
- Electrical height adjustment from 64 - 104 cm
- Lateral up to 40° electric tilting
- Continuously variable V-Top position (possible up to 90°) with one-hand operation by gas spring
- Longitudinal up to 10° electric tilting in both directions possible for (Anti-)Trendelenburg
- Expanded space underneath table top compared to previous version due to optimized design
- Double castors with directional and central locking mechanism by actuating a single pedal
- Detachable and easy-to-clean drain gutter and drip tray for collecting excess liquids

Basic equipment:

OP-System with castors
OP rails (4 pieces on longer side)
Removable laparoscopy cut-outs (2 pieces)
Hand switch
Optional: Joystick



OP-System X-Top

with removable laparoscopy cut-outs



Electrical height adjustable (power cord)

Movement	Item code
movable with centrally-lockable double castors	652.551

Battery-powered height adjustable

Movement	Item code
movable with centrally-lockable double castors	672.551



Power supply: 230V 50/60Hz 3A

Lifting capacity: 120kg