



PURE Veterinary Surgery and Treatment Table

- High-end surgery and treatment table with versatile use in veterinary medicine
- Frame of maximum stability made of first-class Austrian steel; powder coated according to RAL 9007
- Exceptionally hygienic housing of sealed, high-quality stainless steel as well as protective lip gaskets
- Electrical height adjustment by footswitch
- Minimum height of 11 cm resp. 18 cm (incl. stand)
- Maximum height of 112 cm resp. 119 cm (incl. stand)
- Plate (stainless steel or synthetic top) with edge and drain (130 x 60 cm) as well as 14 tie-down hooks
- Available either free-standing with foot or wall-mounted for maximum solidity
- Modular build for retrofittability and versatility
- Top optionally fixed or tiltable (continuously adjustable tilt angle of the top (up to 12°) with comfortable one-hand operation supported by gas spring)
- Timeless elegant and durable design
- Co-developed with world-renowned experts in veterinary medicine and veterinarians
- Lifting capacity up to 120kg with uniform load distribution

**Partially shown OR side rails not included in standard scope of delivery*



Special version with integrated scale

- 0-100 kg in 50 g steps
- Animal weighing function (averaging method)
- Battery-operated (220 hours working hours)
- Automatic shutdown after 3 minutes
- Delivery without batteries (6 pieces AA)



PURE

free-standing with foot



Electrical height adjustable (power cord)

Top	Item Code
fixed stainless steel top	655. 600
tilting stainless steel top	655. 602
fixed synthetic top	655. 610
tilting synthetic top	655. 612
fixed stainless steel top with integrated scale	655. 603
fixed synthetic top with integrated scale	655. 613

PURE

wall-mounted



Electrical height adjustable (power cord)

Top	Item Code
fixed stainless steel top	655. 600-WM
tilting stainless steel top	655. 602-WM
fixed synthetic top	655. 610-WM
tilting synthetic top	655. 612-WM
fixed stainless steel top with integrated scale	655. 603-WM
fixed synthetic top with integrated scale	655. 613-WM

Technical Dimensions

- Dimensions base frame LxBxH: 41 cm (resp. 16 cm without arm) x 55 cm x 145 cm
- Dimensions stand LxBxH: 105 cm x 54 cm (no cap) x 7 cm
- Dimensions top set fixed/tilting LxBxH: 132 cm x 60 cm (64 cm incl. OP rails) x 10 cm
- Dimensions top set with scale LxBxH: 130 cm x 60 cm x 15 cm
- Total dimensions PURE (with stand) LxBxH: 149 cm x 60 cm (64 cm incl. OP rails) x 152 cm
149 cm x 64 cm (with integrated scale) x 152 cm
- Total dimensions PURE (wall-m.) LxBxH: 147 cm x 60 cm (64 cm incl. OP rails) x 145 cm
147 cm x 64 cm (with integrated scale) x 145 cm
- Lifting capacity: 120kg with uniform load distribution
- Height adjustment range incl. stand: from 18 cm to 119 cm
- Height adjustment range (wall-mounted): from 11 cm to 112 cm



Power supply: 230V 50/60Hz 1,5A

Lifting capacity: 120 kg