

Veterinary Dental Series



Better Care For Oral Health



FC Series

FC Series VET dental workstation can be used as a comprehensive dental working platform for tooth polishing, light curing, restoration, cutting and root canal treatment. It can also be used for animal tooth beauty, tooth removal, oral tumor cutting etc., bringing more help to animal welfare, It can also improve the veterinarian's handling of animal oral problems more efficiently.

> Multifunctional dental equipment



Design for animal

Comprehensive & professional

FC Series

Technical Parameters

Specification	FC6-Vet	FC5-Vet	FC4-Vet
Compressor	Oil-Free	Oil-Free	Oil-Free
Power	2.6-2.9A;220/230v; 50/60HZ; 110V(Optional)	2.6-2.9A;220/230v; 50/60HZ; 110V(Optional)	2.6-2.9A;220/230v; 50/60HZ; 110V(Optional)
Out Put Pressure	0.45~0.5MPa	0.45~0.5MPa	0.45~0.5MPa
Gas tank	8L	8L	8L
Sound Level	Less than 60dB	Less than 60dB	Less than 60dB
High Speed Handpiece	300000-340000rpm; Four-point spray; LED light	300000-340000rpm; Four-point spray; LED light	300000-340000rpm; Four-point spray; LED light
Low Speed Handpiece	External water	External water	External water
Low Speed	External water supply	External water supply	External water supply
Motor	Suitable for angle and straight external handpieces	Suitable for angle and straight external handpieces	Suitable for angle and straight external handpieces
Ultrasonic Scaler	With white LED light 28kHz±3kHz 5 scaler tips	With white LED light 28kHz±3kHz 5 scaler tips	With white LED light 28kHz±3kHz 5 scaler tips
Suction handpiece	High volume suction system	High volume suction system	--
LED Light curing unit	Wireless; Wavelength:420-480nm; Depth of cure:3mm/10s;	--	--
Dimension	L*W*H:45 cm * 52 cm * 92 cm		
Weight	About 60kg (Net Weight 45kg)		



Better care for animal 🦷

Product info.

Parameter	FC6-Vet	FC5-Vet	FC4-Vet
Compressor	○	○	○
High Speed Handpiece1	○	○	○
High Speed Handpiece2	optional	--	--
Low Speed Handpiece	○	○	○
Scaler Handpiece	○	○	○
3 Way Syringe	○	○	○
Suction Handpiece	○	○	--
LED Light Curing Unit	○	--	--

(○:YES; --:NO)



Enjoy Health • Enjoy life

www.ylvet.com | Youlan Medical

Shenzhen Youlan Medical Technology Co.,LTD.

Address: 516, Building 4, Baoyunda Logistics Center, Fuhua Community,
Xixiang Street, Bao'an District, Shenzhen, Guangdong Province, China

Tel/WhatsApp: **+0086 188-2019-3720**

Email: sales@ylvet.com



 Facebook



 LinkedIn



 Instagram